

CAVEX

L



Cavex laboratory materials are the natural extension of the dental practice and must therefore meet stringent quality requirements. Thanks to Cavex's close collaboration with dental technicians and prosthetists these materials are also closely matched to the specific requirements of the dental laboratory. For dental use only.

The modelling wax you can read and write with.

Creating a prosthetic model is an art. It's all down to the talent and skills of the dental technician, his or her close collaboration with the dentist and the superior qualities of the modelling wax. For the modelling wax the dental technician accordingly lays down demanding requirements. And Cavex Set Up Modelling Wax really does act like wax in the technician's hands.



Cavex Holland BV
P.O. Box 852
2003 RW Haarlem
The Netherlands



CAVEX

Cavex Holland BV
Harmenjansweg 19
P.O. Box 852
2003 RW Haarlem
The Netherlands

T +31 23 530 77 00
F +31 23 535 64 82
dental@cavex.nl
www.cavex.nl

Predictable behaviour from Cavex Set Up Soft, Regular and Hard Modelling Wax.

- **for full and partial prostheses**
- **low shrinkage stress dimensionally stable**
- **plastic and elastic**
- **modelling and boxing wax**

The secret of a good modelling wax lies in the recipe. Cavex carefully combines paraffins, microwaxes, natural wax and pigments to ensure that the melting points are never too far apart, thus ensuring an even and predictable melting and curing. Because circumstances and personal preferences vary, we offer you a choice of Soft, Regular and Hard. The procedure is the same, but the melting point and the final hardness differ.


The wax plates are of ample size and are usually sufficient for two prostheses. Besides its use as a bite wafer Cavex Set Up Modelling Wax can also be used as a boxing wax because of its excellent adhesion to the impression tray. Thanks to its ideal consistency this wax is popular with dentists for bite registration.

A prosthesis must not only look good, it must be a perfect fit as well.

Aesthetic and technical perfection should go hand-in-hand. To be able to model a wax beautifully and yet still be able to adjust the elements, the wax must remain plastic in the dental technician's hands for a sufficient length of time. Furthermore it retains its dimensional stability. The wax's high degree of elasticity makes it extremely suitable for partial prostheses.

Cavex Set Up Modelling Wax gives an attractive, smooth finish without droplet formation. It offers both an extremely natural colouring - something that both the dentist and the patient will appreciate when it comes to the fitting - and excellent shape retention in the mouth. After curing, the elements will be firmly fixed. The end result is exactly the same as the model on which it is based.

Cavex Set Up Soft, Cavex Set Up Regular and Cavex Set Up Hard are available from your supplier in packs of 500 gram (23 plates) and 2.5 kg (115 plates). For further information on product composition and applications, please contact Cavex, +31 23 530 77 00 - Dental Customer Service or www.cavex.nl



CAVEX FOR DENTAL USE ONLY

ZYW001E01