

PD Fibrapost

Glass fibre composite posts

Properties

Braided glass fibers in a multiaxial arrangement & in a resin bonding matrix (Epoxy resin with approximately 65% glass fiber content) providing high resistance to bending & torsion forces.

Benefits

- Physico-chemical compatibility & high bonding between post, cement & dentine
- Prevents root fracture, oxidation & corrosion
- Flexibility & high resistance
- Aesthetic thanks to its white color
- Color coded for easy identification

Indication

Core build-ups, especially for anterior teeth where aesthetic is of great importance. To be used with corresponding PD spiral bladed tapered reamers with adjustable depth stopper

Packaging

Assortment of:

- 6 posts each N°1 to 4 = 24 posts
- 1 reamer each N°1 to 4 = 4 reamers
- 1 screw driver



FIBRAPOST

