

Oro Clean Chemie AG  
Allmendstrasse 21  
CH-8320 Fehraltorf, Switzerland  
Tel: ++41 1 954 80 30  
Fax: ++41 1 954 80 40

# Surface disinfection

## Quick disinfection

|                                   |  |
|-----------------------------------|--|
| <b>PRODUCT:</b>                   | <b>DENTIRO<sup>®</sup> Sensitive</b><br>Special product for the quick disinfection of medical equipment, unit surfaces and instruments (in the sense of commission directive 93/42/EEC). Can be used for other small surfaces. Aldehyde free. Easily biodegradable. Ready-to-use solution.   |
| <b>Use:</b>                       | Moisten the surfaces and the objects completely, allow reacting and dry slightly if necessary. Use only for alcohol compatible materials. With sensitive materials test compatibility before use.  |
| <b>Activities:</b>                | <input type="checkbox"/> Bactericide, fungicide within 30 seconds.<br><input type="checkbox"/> Tuberculocide within 30 seconds.<br><input type="checkbox"/> Hepatitis B <sup>1)</sup> incl. HIV within 30 seconds.<br><input type="checkbox"/> Hepatitis C <sup>2)</sup> within 30 seconds.<br>1) HBsAg – antigen test; 2) BVDV – surrogate virus test   |
| <b>Composition in 100 g:</b>      | 20 g ethanol, 28 g 1-propanol, 0.056 g quaternary ammonium compounds.  |
| <b>Physical properties:</b>       | Clear, colourless liquid with alcoholic odour.<br>Flash point: 26 °C (DIN 51 755)<br>Density (20 °C): ~ 0.906 g/ml (g/cm <sup>3</sup> )  |
| <b>Toxicological information:</b> | LD <sub>50</sub> - Value: > 5000 mg/kg   |
| <b>Registrations:</b>             | CE 0483<br>Germany: DGHM listed  |
| <b>Shelf life:</b>                | 36 months (Lot and expiry-date on containment)   |
| <b>Precautionary measures:</b>    | <b>R10</b> Flammable. <b>R41</b> Risk of serious damage to eyes. <b>S2</b> Keep out of reach of children. <b>S16</b> Keep away from sources of ignition. No smoking. <b>S23</b> Do not breathe spray. <b>S25</b> Avoid contact with eyes. <b>S25</b> In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. <b>S51</b> Use only in well-ventilated areas. |
| <b>Test reports:</b>              | Prof. Dr. med. R. Schubert, Frankfurt, D, 2001<br>Prof. Dr. med. H.-P. Werner, Schwerin, D, 2002, 2001<br>Dr. J. Steinmann, Bremen, D, 2001  |

These information contained herein are believed to be correct and correspond to the best knowledge available. However, Oro Clean Chemie AG, Allmendstrasse 21 CH-8320 Fehraltorf, Switzerland cannot take over the responsibility regarding its accuracy or the completeness of these information.