




Impression materials are the starting point of many dental operations. High-quality impressions give the dentist confidence in a good end result. Cavex offers a well-tried and reliable range of alginates, hydrocolloid range of alginates, hydrocolloid and auxiliary materials developed according to the requirements of the dental practice. For dental use only.

Cavex CA37, because your first impression must also be the best.

Alginate is the most versatile impression material available in the dental practice. With Cavex CA37 you can be sure of obtaining the accuracy that can be decisive for a successful treatment. Tested in practice for 50 years, Cavex CA37 can boast tens of millions of perfect impressions worldwide per year resulting in just as many successful treatments. Is this a good enough reference for your practice?





CAVEX

Cavex Holland BV
Harmerjansweg 19
P.O. Box 852
2003 RW Haarlem
The Netherlands

T +31 23 530 77 00
F +31 23 535 64 82
dental@cavex.nl
www.cavex.nl

Better to get an accurate impression at once than have to correct it later.

- **ISO 1563 / ADA 18**

Whether it's a first impression, a functional impression, an impression for a countermodel, a study model...what you want is a highly detailed impression that will supply an accurate model of the situation in the mouth. After all, any tolerance will wreak its revenge in later stages of the treatment.

- **worldwide references**

This quality of impression is guaranteed by Cavex CA37 and confirmed by tests according to the international institutes ISO and ADA, who have laid down the highest standards for detail reproduction and readability. The high degree of compatibility with plaster gives a smooth surface and ensures a perfect, true-to-life model.

- **dust-free dosing**

- **elastic, reversible, tear-resistant**

Fluid under pressure, stable on the tray, the ideal working procedure.

Getting the dosing and mixing right helps ensure a perfect result. So it's good that Cavex CA37 is dust-free and easy to mix. The rapid moisture absorption encourages a full reaction with water and eliminates powder clumps.

- **short curing time**

- **compatible with plaster, smooth surface**

The alginate mass has the ideal consistency - it's fluid under pressure but stable on the tray. After curing, the impression remains elastic and is easily removed from the mouth despite undercuts. The impression immediately springs back into its original shape.

The patient will appreciate the short curing time and the subtle peppermint flavour. The total processing time, including mixing, is just 3,30 minutes for Normal Set and 2,30 minutes for Fast Set.

Cavex CA37 Normal Set (Superior Pink) and Fast Set are available from your local dental supplier. For further information about product composition and applications, please contact Cavex, +31 23 530 77 00 - Dental Customer Service.



CAVEX FOR DENTAL USE ONLY

ZYEC2000/2-0