

For the Banding Technique



Ortho-Cast™ · single round

Round tubes 0,9 mm / 36 · non-convertible · with vertical slot · for lip bumper / face bow technique

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
17-16/36-37	-	6°	0,3 mm	•	4,0 mm	724-080-00	10 pieces
26-27/47-46	-	6°	0,3 mm	•	4,0 mm	724-081-00	10 pieces

CE 0483



Ortho-Cast™ 18 · single rectangular

Slot: 0,46 x 0,64 mm / 18 x 25 · non-convertible

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
17-16	- 10°	8°	0,3 mm	•	4,0 mm	724-016-00	10 pieces
26-27	- 10°	8°	0,3 mm	•	4,0 mm	724-017-00	10 pieces

CE 0483



Ortho-Cast™ 18 · single rectangular

Slot: 0,46 x 0,64 mm / 18 x 25 · convertible

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
17-16	- 10°	8°	0,3 mm	•	4,0 mm	724-048-00	10 pieces
26-27	- 10°	8°	0,3 mm	•	4,0 mm	724-049-00	10 pieces

CE 0483



Ortho-Cast™ 18 · single rectangular

Slot: 0,46 x 0,64 mm / 18 x 25 · non-convertible

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
47-46	- 25°	8°	0,3 mm	•	4,0 mm	724-020-00	10 pieces
36-37	- 25°	8°	0,3 mm	•	4,0 mm	724-021-00	10 pieces
47-46	- 30°	12°	0,3 mm	•	4,0 mm	724-014-00	10 pieces
36-37	- 30°	12°	0,3 mm	•	4,0 mm	724-015-00	10 pieces

CE 0483



Ortho-Cast™ 18 · single rectangular

Slot: 0,46 x 0,64 mm / 18 x 25 · convertible

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
47-46	- 25°	8°	0,3 mm	•	4,0 mm	724-040-00	10 pieces
36-37	- 25°	8°	0,3 mm	•	4,0 mm	724-041-00	10 pieces

CE 0483



Ortho-Cast™ 18 · single rectangular

Slot: 0,46 x 0,64 mm / 18 x 25 · non-convertible

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
17-16/36-37	0°	6°	0,3 mm	•	4,0 mm	724-092-00	10 pieces
26-27/47-46	0°	6°	0,3 mm	•	4,0 mm	724-093-00	10 pieces
17-16/26-27/47-46/36-37	0°	6°	0,3 mm	•	4,0 mm	724-090-00	10 pieces

CE 0483



Ortho-Cast™ 18 · single rectangular

Slot: 0,46 x 0,64 mm / 18 x 25 · convertible

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
17-16/47-46	0°	8°	0,3 mm	•	4,0 mm	724-044-00	10 pieces
26-27/36-37	0°	8°	0,3 mm	•	4,0 mm	724-045-00	10 pieces

CE 0483

Ortho-Cast™ 18 · single rectangular

Slot: 0,46 x 0,64 mm / 18 x 25 · non-convertible · distal extension

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
17-16/36-37	0°	0°	0,3 mm	•	6,0 mm	724-082-00	10 pieces
26-27/47-46	0°	0°	0,3 mm	•	6,0 mm	724-083-00	10 pieces
17-16/36-37	0°	6°	0,3 mm	•	6,0 mm	724-086-00	10 pieces
26-27/47-46	0°	6°	0,3 mm	•	6,0 mm	724-087-00	10 pieces



CE 0483

Ortho-Cast™ 18 · double rectangular

Slot: 0,46 x 0,64 mm / 18 x 25 · non-convertible

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
17-16	- 9°	10°	0,3 mm	•	4,0 mm	724-066-00	10 pieces
26-27	- 9°	10°	0,3 mm	•	4,0 mm	724-067-00	10 pieces



CE 0483

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
47-46	- 22°	12°	0,3 mm	•	4,0 mm	724-064-00	10 pieces
36-37	- 22°	12°	0,3 mm	•	4,0 mm	724-065-00	10 pieces



CE 0483

Ortho-Cast™ 18 · double rectangular

Slot: 0,46 x 0,64 mm / 18 x 25 · non-convertible · gingival auxiliary tube

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
47-46	- 30°	4°	0,3 mm	•	3,5 mm	724-208-00	10 pieces
36-37	- 30°	4°	0,3 mm	•	3,5 mm	724-209-00	10 pieces



CE 0483

Ortho-Cast™ 18 · double rectangular

Slot: 0,46 x 0,64 mm / 18 x 25 · non-convertible · with vertical hook

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
47-46/36-37	- 30°	0°	0,3 mm	•	4,0 mm	724-062-00	10 pieces



CE 0483

Ortho-Cast™ 18 · double rectangular

Slot: 0,46 x 0,64 mm / 18 x 25 · convertible

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
47-46	0°	0°	0,3 mm	•	4,0 mm	724-002-00	10 pieces
36-37	0°	0°	0,3 mm	•	4,0 mm	724-003-00	10 pieces
47-46	- 25°	8°	0,3 mm	•	4,0 mm	724-056-00	10 pieces
36-37	- 25°	8°	0,3 mm	•	4,0 mm	724-057-00	10 pieces
47-46	- 25°	12°	0,3 mm	•	4,0 mm	724-012-00	10 pieces
36-37	- 25°	12°	0,3 mm	•	4,0 mm	724-013-00	10 pieces



CE 0483

Ortho-Cast™ 18 · double rectangular

Slot: 0,46 x 0,64 mm / 18 x 25 · non-convertible · with universal vertical hook

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
17-16/26-27/47-46/36-37	0°	0°	0,3 mm	•	3,5 mm	724-018-00	10 pieces



CE 0483

For the Banding Technique



Ortho-Cast™ 18 · double combination

Slot: 0,46 x 0,64 mm / 18 x 25 · tube 1,15 mm / 45 (occlusal) · non-convertible · for lip bumper / face bow technique

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
17-16	- 14°	14°	0,3 mm	•	4,0 mm	724-204-00	10 pieces
26-27	- 14°	14°	0,3 mm	•	4,0 mm	724-205-00	10 pieces

CE 0483



Ortho-Cast™ 18 · double combination

Slot: 0,46 x 0,64 mm / 18 x 25 · tube 1,15 mm / 45 (occlusal) · convertible · for lip bumper / face bow technique

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
17-16	0°	8°	0,3 mm	•	4,0 mm	724-004-00	10 pieces
26-27	0°	8°	0,3 mm	•	4,0 mm	724-005-00	10 pieces
17-16	- 10°	8°	0,3 mm	•	4,0 mm	724-028-00	10 pieces
26-27	- 10°	8°	0,3 mm	•	4,0 mm	724-029-00	10 pieces
17-16	- 10°	15°	0,3 mm	•	4,0 mm	724-058-00	10 pieces
26-27	- 10°	15°	0,3 mm	•	4,0 mm	724-059-00	10 pieces

CE 0483



Ortho-Cast™ 18 · double combination

Slot: 0,46 x 0,64 mm / 18 x 25 · tube 1,15 mm / 45 (gingival) · convertible · for lip bumper / face bow technique

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
47-46	- 25°	8°	0,3 mm	•	4,0 mm	724-036-00	10 pieces
36-37	- 25°	8°	0,3 mm	•	4,0 mm	724-037-00	10 pieces

CE 0483



Ortho-Cast™ 18 · triple combination

Slot: 0,46 x 0,64 mm / 18 x 25 · round tube 1,15 mm / 45 (occlusal) · non-convertible · for lip bumper / face bow technique

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
17-16	- 14°	14°	0,3 mm	•	4,0 mm	724-200-00	10 pieces
26-27	- 14°	14°	0,3 mm	•	4,0 mm	724-201-00	10 pieces

CE 0483



Ortho-Cast™ 18 · triple combination

Slot: 0,46 x 0,64 mm / 18 x 25 · round tube 1,15 mm / 45 (central) · non-convertible · for lip bumper / face bow technique

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
17-16	- 9°	10°	0,3 mm	•	4,0 mm	724-068-00	10 pieces
26-27	- 9°	10°	0,3 mm	•	4,0 mm	724-069-00	10 pieces

CE 0483



Ortho-Cast™ 18 · triple combination

Slot: 0,46 x 0,64 mm / 18 x 25 · round tube 1,15 mm / 45 (occlusal) · convertible · for lip bumper / face bow technique

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
17-16	0°	8°	0,3 mm	•	4,0 mm	724-006-00	10 pieces
26-27	0°	8°	0,3 mm	•	4,0 mm	724-007-00	10 pieces
17-16	0°	15°	0,3 mm	•	4,0 mm	724-008-00	10 pieces
26-27	0°	15°	0,3 mm	•	4,0 mm	724-009-00	10 pieces
17-16	- 10°	15°	0,3 mm	•	4,0 mm	724-010-00	10 pieces
26-27	- 10°	15°	0,3 mm	•	4,0 mm	724-011-00	10 pieces

CE 0483



Ortho-Cast™ 18 · triple combination

Slot: 0,46 x 0,64 mm / 18 x 25 · round tube 1,15 mm / 45 (central) · non-convertible · for lip bumper / face bow technique

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
47-46	- 25°	8°	0,3 mm	•	4,0 mm	724-060-00	10 pieces
36-37	- 25°	8°	0,3 mm	•	4,0 mm	724-061-00	10 pieces

CE 0483

Ortho-Cast™ 22 · single rectangular

Slot: 0,56 x 0,70 mm / 22 x 28 · non-convertible

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
17-16	- 10°	8°	0,3 mm	•	4,0 mm	724-116-00	10 pieces
26-27	- 10°	8°	0,3 mm	•	4,0 mm	724-117-00	10 pieces
17-16	- 14°	10°	0,3 mm	•	4,0 mm	724-214-00	10 pieces
26-27	- 14°	10°	0,3 mm	•	4,0 mm	724-215-00	10 pieces



CE 0483

Ortho-Cast™ 22 · single rectangular

Slot: 0,56 x 0,70 mm / 22 x 28 · convertible

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
17-16	- 10°	8°	0,3 mm	•	4,0 mm	724-148-00	10 pieces
26-27	- 10°	8°	0,3 mm	•	4,0 mm	724-149-00	10 pieces



CE 0483

Ortho-Cast™ 22 · single rectangular

Slot: 0,56 x 0,70 mm / 22 x 28 · non-convertible

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
47-46	- 10°	0°	0,3 mm	•	4,0 mm	724-218-00	10 pieces
36-37	- 10°	0°	0,3 mm	•	4,0 mm	724-219-00	10 pieces
47-46	- 25°	8°	0,3 mm	•	4,0 mm	724-120-00	10 pieces
36-37	- 25°	8°	0,3 mm	•	4,0 mm	724-121-00	10 pieces
47-46	- 30°	12°	0,3 mm	•	4,0 mm	724-114-00	10 pieces
36-37	- 30°	12°	0,3 mm	•	4,0 mm	724-115-00	10 pieces



CE 0483

Ortho-Cast™ 22 · single rectangular

Slot: 0,56 x 0,70 mm / 22 x 28 · convertible

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
47-46	- 20°	0°	0,3 mm	•	4,0 mm	724-228-00	10 pieces
36-37	- 20°	0°	0,3 mm	•	4,0 mm	724-229-00	10 pieces
47-46	- 25°	8°	0,3 mm	•	4,0 mm	724-140-00	10 pieces
36-37	- 25°	8°	0,3 mm	•	4,0 mm	724-141-00	10 pieces



CE 0483

Ortho-Cast™ 22 · single rectangular

Slot: 0,56 x 0,70 mm / 22 x 28 · non-convertible

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
17-16/36-37	0°	6°	0,3 mm	•	4,0 mm	724-094-00	10 pieces
26-27/47-46	0°	6°	0,3 mm	•	4,0 mm	724-095-00	10 pieces
17-16/26-27/47-46/36-37	0°	6°	0,3 mm	•	4,0 mm	724-091-00	10 pieces



CE 0483

Ortho-Cast™ 22 · single rectangular

Slot: 0,56 x 0,70 mm / 22 x 28 · convertible

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
17-16/36-37	0°	8°	0,3 mm	•	4,0 mm	724-144-00	10 pieces
26-27/47-46	0°	8°	0,3 mm	•	4,0 mm	724-145-00	10 pieces



CE 0483

For the Banding Technique



Ortho-Cast™ 22 · single rectangular

Slot: 0,56 x 0,70 mm / 22 x 28 · non-convertible · distal extension

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
17-16/36-37	0°	0°	0,3 mm	•	6,0 mm	724-084-00	10 pieces
26-27/47-46	0°	0°	0,3 mm	•	6,0 mm	724-085-00	10 pieces
17-16/36-37	0°	6°	0,3 mm	•	6,0 mm	724-088-00	10 pieces
26-27/47-46	0°	6°	0,3 mm	•	6,0 mm	724-089-00	10 pieces

CE 0483



Ortho-Cast™ 22 · double rectangular

Slot: 0,56 x 0,70 mm / 22 x 28 · non-convertible

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
17-16	- 9°	10°	0,3 mm	•	4,0 mm	724-166-00	10 pieces
26-27	- 9°	10°	0,3 mm	•	4,0 mm	724-167-00	10 pieces

CE 0483



Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
47-46	- 22°	12°	0,3 mm	•	4,0 mm	724-164-00	10 pieces
36-37	- 22°	12°	0,3 mm	•	4,0 mm	724-165-00	10 pieces

CE 0483



Ortho-Cast™ 22 · double rectangular

Slot: 0,56 x 0,70 mm / 22 x 28 · non-convertible · gingival auxiliary tube

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
47-46	- 30°	4°	0,3 mm	•	3,5 mm	724-210-00	10 pieces
36-37	- 30°	4°	0,3 mm	•	3,5 mm	724-211-00	10 pieces

CE 0483



Ortho-Cast™ 22 · double rectangular

Slot: 0,56 x 0,70 mm / 22 x 28 · non-convertible · with vertical hook

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
47-46/36-37	- 30°	0°	0,3 mm	•	4,0 mm	724-063-00	10 pieces

CE 0483



Ortho-Cast™ 22 · double rectangular

Slot: 0,56 x 0,70 mm / 22 x 28 · convertible

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
47-46	0°	0°	0,3 mm	•	4,0 mm	724-102-00	10 pieces
36-37	0°	0°	0,3 mm	•	4,0 mm	724-103-00	10 pieces
47-46	- 25°	8°	0,3 mm	•	4,0 mm	724-156-00	10 pieces
36-37	- 25°	8°	0,3 mm	•	4,0 mm	724-157-00	10 pieces
47-46	- 25°	12°	0,3 mm	•	4,0 mm	724-112-00	10 pieces
36-37	- 25°	12°	0,3 mm	•	4,0 mm	724-113-00	10 pieces

CE 0483



Ortho-Cast™ 22 · double rectangular

Slot: 0,56 x 0,70 mm / 22 x 28 · non-convertible · with universal vertical hook

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
17-16/26-27/47-46/36-37	0°	0°	0,3 mm	•	3,5 mm	724-019-00	10 pieces

CE 0483

Ortho-Cast™ 22 · double combination

Slot: 0,56 x 0,70 mm / 22 x 28 · tube 1,15 mm / 45 (occlusal) · convertible · for lip bumper / face bow technique

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
17-16	0°	8°	0,3 mm	•	4,0 mm	724-104-00	10 pieces
26-27	0°	8°	0,3 mm	•	4,0 mm	724-105-00	10 pieces
17-16	- 10°	8°	0,3 mm	•	4,0 mm	724-030-00	10 pieces
26-27	- 10°	8°	0,3 mm	•	4,0 mm	724-031-00	10 pieces
17-16	- 10°	15°	0,3 mm	•	4,0 mm	724-158-00	10 pieces
26-27	- 10°	15°	0,3 mm	•	4,0 mm	724-159-00	10 pieces



CE 0483

Ortho-Cast™ 22 · double combination

Slot: 0,56 x 0,70 mm / 22 x 28 · tube 1,15 mm / 45 (occlusal) · non-convertible · for lip bumper / face bow technique

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
17-16	- 14°	14°	0,3 mm	•	4,0 mm	724-206-00	10 pieces
26-27	- 14°	14°	0,3 mm	•	4,0 mm	724-207-00	10 pieces



CE 0483

Ortho-Cast™ 22 · double combination

Slot: 0,56 x 0,70 mm / 22 x 28 · tube 1,15 mm / 45 (gingival) · convertible · for lip bumper / face bow technique

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
17-16	- 10°	8°	0,3 mm	•	4,0 mm	724-150-00	10 pieces
26-27	- 10°	8°	0,3 mm	•	4,0 mm	724-151-00	10 pieces



CE 0483

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
47-46	- 25°	8°	0,3 mm	•	4,0 mm	724-038-00	10 pieces
36-37	- 25°	8°	0,3 mm	•	4,0 mm	724-039-00	10 pieces



CE 0483

Ortho-Cast™ 22 · triple combination

Slot: 0,56 x 0,70 mm / 22 x 28 · round tube 1,15 mm / 45 (occlusal) · non-convertible · for lip bumper / face bow technique

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
17-16	- 14°	14°	0,3 mm	•	4,0 mm	724-202-00	10 pieces
26-27	- 14°	14°	0,3 mm	•	4,0 mm	724-203-00	10 pieces



CE 0483

Ortho-Cast™ 22 · triple combination

Slot: 0,56 x 0,70 mm / 22 x 28 · round tube 1,15 mm / 45 (central) · non-convertible · for lip bumper / face bow technique

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
17-16	- 9°	10°	0,3 mm	•	4,0 mm	724-070-00	10 pieces
26-27	- 9°	10°	0,3 mm	•	4,0 mm	724-071-00	10 pieces



CE 0483

Ortho-Cast™ 22 · triple combination

Slot: 0,56 x 0,70 mm / 22 x 28 · round tube 1,15 mm / 45 (occlusal) · convertible · for lip bumper / face bow technique

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
17-16	0°	8°	0,3 mm	•	4,0 mm	724-106-00	10 pieces
26-27	0°	8°	0,3 mm	•	4,0 mm	724-107-00	10 pieces
17-16	- 10°	8°	0,3 mm	•	4,0 mm	724-162-00	10 pieces
26-27	- 10°	8°	0,3 mm	•	4,0 mm	724-163-00	10 pieces
17-16	- 10°	15°	0,3 mm	•	4,0 mm	724-110-00	10 pieces
26-27	- 10°	15°	0,3 mm	•	4,0 mm	724-111-00	10 pieces
17-16	- 14°	10°	0,3 mm	•	4,0 mm	724-226-00	10 pieces
26-27	- 14°	10°	0,3 mm	•	4,0 mm	724-227-00	10 pieces



CE 0483

For the Banding Technique



Ortho-Cast™ 22 · triple combination

Slot: 0,56 x 0,70 mm / 22 x 28 · round tube 1,15 mm / 45 (central) · non-convertible · for lip bumper / face bow technique

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
47-46	- 25°	8°	0,3 mm	•	4,0 mm	724-160-00	10 pieces
36-37	- 25°	8°	0,3 mm	•	4,0 mm	724-161-00	10 pieces

CE 0483



Ortho · single round

round tube 0,9 mm / 36 · non-convertible · for Light Wire Technique · with vertical slot · for lip bumper / face bow technique

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
17-16/36-37	-	6°	0,3 mm	•	6,0 mm	725-033-00	10 pieces
26-27/47-46	-	6°	0,3 mm	•	6,0 mm	725-034-00	10 pieces

CE 0483



Ortho · single round

round tube 1,15 mm / 45 · non-convertible · for lip bumper / face bow technique

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
17-16/26-27/47-46/36-37	-	0°	0,3 mm		6,0 mm	725-206-00	10 pieces

CE 0483



Ortho · double round

round tube 0,9 mm / 36 · round tube 1,15 mm / 45 · distal extension and vertical slot · for lip bumper / face bow technique

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
17-16/36-37	-	6°	0,3 mm	•	4,6 mm	726-009-00	10 pieces
26-27/47-46	-	6°	0,3 mm	•	4,6 mm	726-010-00	10 pieces

CE 0483



Ortho · double round

round tube 0,9 mm / 36 · round tube 1,15 mm / 45 · for lip bumper / face bow technique

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
17-16/36-37	-	0°	0,3 mm	•	4,6 mm	726-012-00	10 pieces
26-27/47-46	-	0°	0,3 mm	•	4,6 mm	726-013-00	10 pieces
17-16/26-27/47-46/36-37	-	0°	0,3 mm		4,6 mm	726-011-00	10 pieces

CE 0483



Ortho 18 · single rectangular

Slot: 0,46 x 0,64 mm / 18 x 25 · non-convertible

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
17-16/26-27/47-46/36-37	0°	6°	0,3 mm		4,0 mm	725-310-00	10 pieces

CE 0483



Ortho 18 · single rectangular

Slot: 0,46 x 0,64 mm / 18 x 25 · non-convertible · distal extension

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
17-16/36-37	0°	0°	0,3 mm	•	6,0 mm	725-164-00	10 pieces
26-27/47-46	0°	0°	0,3 mm	•	6,0 mm	725-165-00	10 pieces
17-16/36-37	0°	6°	0,3 mm	•	6,0 mm	725-168-00	10 pieces
26-27/47-46	0°	6°	0,3 mm	•	6,0 mm	725-169-00	10 pieces

CE 0483



Ortho 18 · single rectangular

Slot: 0,46 x 0,64 mm / 18 x 25 · non-convertible

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
17-16/36-37	0°	6°	0,3 mm	•	4,0 mm	725-308-00	10 pieces
26-27/47-46	0°	6°	0,3 mm	•	4,0 mm	725-309-00	10 pieces

CE 0483

Ortho 18 · double combination

Slot: 0,46 x 0,64 mm / 18 x 25 · non-convertible · for lip bumper / face bow technique

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
17-16	0°	0°	0,3 mm	•	4,6 mm	727-365-00	10 pieces
26-27	0°	0°	0,3 mm	•	4,6 mm	727-366-00	10 pieces
17-16	0°	6°	0,3 mm	•	4,6 mm	727-367-00	10 pieces
26-27	0°	6°	0,3 mm	•	4,6 mm	727-368-00	10 pieces
17-16/26-27	0°	0°	0,3 mm		4,6 mm	727-265-00	10 pieces



CE 0483

Ortho 18 · double combination

Slot: 0,46 x 0,64 mm / 18 x 25 · round tube 1,15 mm / 45 (gingival) · non-convertible · for lip bumper / face bow technique

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
47-46	0°	0°	0,3 mm	•	4,6 mm	727-267-00	10 pieces
36-37	0°	0°	0,3 mm	•	4,6 mm	727-268-00	10 pieces



CE 0483

Ortho 18 · triple combination

Slot: 0,46 x 0,64 mm / 18 x 25 · round tube 1,15 mm / 45 (gingival) · non-convertible · for lip bumper / face bow technique

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
17-16	0°	0°	0,3 mm	•	4,6 mm	726-208-00	10 pieces
26-27	0°	0°	0,3 mm	•	4,6 mm	726-209-00	10 pieces



CE 0483

Ortho 18 · triple combination

Slot: 0,46 x 0,64 mm / 18 x 25 · round tube 1,15 mm / 45 (central) · non-convertible · for lip bumper / face bow technique

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
17-16/26-27	0°	0°	0,3 mm		4,6 mm	726-201-00	10 pieces
17-16/26-27	0°	6°	0,3 mm		4,6 mm	726-203-00	10 pieces
17-16	0°	0°	0,3 mm	•	4,6 mm	726-206-00	10 pieces
17-16	0°	6°	0,3 mm	•	4,6 mm	726-204-00	10 pieces
26-27	0°	0°	0,3 mm	•	4,6 mm	726-207-00	10 pieces
26-27	0°	6°	0,3 mm	•	4,6 mm	726-205-00	10 pieces



CE 0483

Ortho 18 · triple combination

Slot: 0,46 x 0,64 mm / 18 x 25 · round tube 1,15 mm / 45 (central) · non-convertible · with mesial/distal hook · for lip bumper / face bow technique

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
17-16/26-27	0°	0°	0,3 mm	•	4,6 mm	726-202-00	10 pieces



CE 0483

Ortho 22 · single rectangular

Slot: 0,56 x 0,70 mm / 22 x 28 · non-convertible · distal extension

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
17-16/36-37	0°	0°	0,3 mm	•	6,0 mm	725-166-00	10 pieces
26-27/47-46	0°	0°	0,3 mm	•	6,0 mm	725-167-00	10 pieces
17-16/36-37	0°	6°	0,3 mm	•	6,0 mm	725-170-00	10 pieces
26-27/47-46	0°	6°	0,3 mm	•	6,0 mm	725-171-00	10 pieces



CE 0483

Ortho 22 · single rectangular

Slot: 0,56 x 0,70 mm / 22 x 28 · non-convertible

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
17-16/36-37	0°	6°	0,3 mm	•	4,0 mm	725-306-00	10 pieces
26-27/47-46	0°	6°	0,3 mm	•	4,0 mm	725-307-00	10 pieces
17-16/26-27/47-46/36-37	0°	6°	0,3 mm		4,0 mm	725-311-00	10 pieces



CE 0483

For the Banding Technique



Ortho 22 · double combination

Slot: 0,56 x 0,70 mm / 22 x 28 · **non-convertible** · for lip bumper / face bow technique

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
17-16	0°	0°	0,3 mm	•	4,6 mm	727-385-00	10 pieces
26-27	0°	0°	0,3 mm	•	4,6 mm	727-386-00	10 pieces
17-16/26-27	0°	0°	0,3 mm		4,6 mm	727-266-00	10 pieces
17-16	0°	6°	0,3 mm	•	4,6 mm	727-387-00	10 pieces
26-27	0°	6°	0,3 mm	•	4,6 mm	727-388-00	10 pieces

CE 0483



Ortho 22 · double combination

Slot: 0,56 x 0,70 mm / 22 x 28 · **round tube 1,15 mm / 45 (occlusal)** · **non-convertible** · **distal extension** · for lip bumper / face bow technique

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
17-16	0°	6°	0,3 mm	•	4,6 mm	727-413-00	10 pieces
26-27	0°	6°	0,3 mm	•	4,6 mm	727-415-00	10 pieces
17-16	0°	0°	0,7 mm	•	4,6 mm	726-210-00	10 pieces
26-27	0°	0°	0,7 mm	•	4,6 mm	726-211-00	10 pieces

CE 0483



Ortho 22 · triple combination

Slot: 0,56 x 0,70 mm / 22 x 28 · **round tube 1,15 mm / 45 (central)** · **non-convertible** · for lip bumper / face bow technique

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
17-16/26-27	0°	0°	0,3 mm		4,6 mm	726-401-00	10 pieces
17-16	0°	6°	0,3 mm	•	4,6 mm	726-404-00	10 pieces
26-27	0°	6°	0,3 mm	•	4,6 mm	726-405-00	10 pieces

CE 0483



Standard · single round

round tube 0,8 mm / 31

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
17-16/26-27/47-46/36-37	-	0°	0,3 mm		6,0 mm	725-002-00	10 pieces

CE 0483



Standard · single round

round tube 0,9 mm / 36

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
17-16/26-27/47-46/36-37	-	0°	0,3 mm		6,0 mm	725-003-00	10 pieces
17-16/26-27/47-46/36-37	-	0°	0,3 mm		4,6 mm	725-023-00	10 pieces

CE 0483



Standard 18 · double rectangular

Slot: 0,46 x 0,64 mm / 18 x 25 · **non-convertible**

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
47-46/36-37	0°	0°	0,3 mm		4,6 mm	726-100-00	10 pieces

CE 0483



Standard 18 · double rectangular

Slot: 0,46 x 0,64 mm / 18 x 25 · **non-convertible** · **with vertical hook**

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
47-46/36-37	0°	0°	0,3 mm	•	4,6 mm	726-101-00	10 pieces

CE 0483

Mollin-Light Wire · single rectangular

Slot: 0,56 x 0,70 mm / 22 x 28 · non-convertible

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
17-16/36-37	-	6°	0,3 mm	•	7,0 mm	725-505-00	10 pieces
26-27/47-46	-	6°	0,3 mm	•	7,0 mm	725-506-00	10 pieces



CE 0483

Mollin-Light Wire · double combination

Slot: 0,46 x 0,64 mm / 18 x 25 · round tube 1,2 mm / 48 · non-convertible · for lip bumper / face bow technique

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
17-16/36-37	-	6°	0,3 mm	•	6,35 mm	725-500-00	10 pieces
26-27/47-46	-	6°	0,3 mm	•	6,35 mm	725-501-00	10 pieces



CE 0483

Preweld - Lingual Attachments



Lingual / Palatal Sheaths

Slot: 2 x 0,9 mm / 2 x 36

For Orthorama®, Lingual / Palatal Bars, Goshgarian, Quad Helix. Open, without hook.

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
17-16/26-27/47-46/36-37	-	0°	0,25 mm		3,5 mm	728-010-00	10 pieces

CE 0483



Lingual / Palatal Sheaths

Slot: 2 x 0,9 mm / 2 x 36

For Orthorama®, Lingual / Palatal Bars, Goshgarian, Quad Helix. Closed, without hook.

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
17-16/26-27/47-46/36-37	-	0°	0,25 mm		3,5 mm	728-110-00	10 pieces

CE 0483



Lingual / Palatal Sheaths

Slot: 2 x 0,9 mm / 2 x 36

For Orthorama®, Lingual / Palatal Bars, Goshgarian, Quad Helix. Open, with hook.

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
17-16/36-37	-	0°	0,25 mm	•	3,5 mm	728-011-00	10 pieces
47-46/26-27	-	0°	0,25 mm	•	3,5 mm	728-012-00	10 pieces

CE 0483



Lingual / Palatal Sheaths

Slot: 2 x 0,9 mm / 2 x 36

For Orthorama®, Lingual / Palatal Bars, Goshgarian, Quad Helix. Closed, with hook.

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
17-16/36-37	-	0°	0,25 mm	•	3,5 mm	728-111-00	10 pieces
47-46/26-27	-	0°	0,25 mm	•	3,5 mm	728-112-00	10 pieces

CE 0483



750-002-00	10 pieces
------------	-----------

CE 0483

Banding Technique Lingual Buttons

For elastics or Dentalastics®.
Round base, long neck.



750-101-00	10 pieces
------------	-----------

CE 0483

Broad base, short neck.



750-201-00	100 pieces
------------	------------

CE 0483

Banding Technique Band Seating Lugs

Suitable for all preformed bands.
On strips

Preweld - Lingual Attachments

Banding Technique Lingual Double Hooks

For additional elastics or tension springs.
Narrow, for cuspid and bicuspid bands.
On strips

750-602-00	10 pieces
------------	-----------

CE 0483



Single

750-602-10	10 pieces
------------	-----------

CE 0483



Banding Technique Lingual Double Hooks

For additional elastics or tension springs.
Wide, for molar bands.
On strips

750-603-00	10 pieces
------------	-----------

CE 0483



Single

750-603-10	10 pieces
------------	-----------

CE 0483



Banding Technique Lingual Button Hooks

For elastics, ligatures, etc.

750-701-00	10 pieces
------------	-----------

CE 0483



Banding Technique Lingual Hooks

For elastics, ligatures, etc.
Single

750-401-00	100 pieces
------------	------------

CE 0483



For the Bonding Technique



Ortho-Cast™ 18 · single rectangular

Slot: 0,46 x 0,64 mm / 18 x 25 · non-convertible

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
17-16	- 10°	8°	0,3 mm	•	4,0 mm	724-016-51	10 pieces
26-27	- 10°	8°	0,3 mm	•	4,0 mm	724-017-51	10 pieces

CE 0483



Ortho-Cast™ 18 · single rectangular

Slot: 0,46 x 0,64 mm / 18 x 25 · convertible

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
17-16	- 10°	8°	0,3 mm	•	4,0 mm	724-048-51	10 pieces
26-27	- 10°	8°	0,3 mm	•	4,0 mm	724-049-51	10 pieces

CE 0483



Ortho-Cast™ 18 · single rectangular

Slot: 0,46 x 0,64 mm / 18 x 25 · non-convertible

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
47-46	- 25°	8°	0,3 mm	•	4,0 mm	724-020-51	10 pieces
36-37	- 25°	8°	0,3 mm	•	4,0 mm	724-021-51	10 pieces
47-46	- 30°	12°	0,3 mm	•	4,0 mm	724-014-51	10 pieces
36-37	- 30°	12°	0,3 mm	•	4,0 mm	724-015-51	10 pieces

CE 0483



Ortho-Cast™ 18 · single rectangular

Slot: 0,46 x 0,64 mm / 18 x 25 · convertible

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
47-46	- 25°	8°	0,3 mm	•	4,0 mm	724-040-51	10 pieces
36-37	- 25°	8°	0,3 mm	•	4,0 mm	724-041-51	10 pieces

CE 0483



Ortho-Cast™ 18 · single rectangular

Slot: 0,46 x 0,64 mm / 18 x 25 · non-convertible

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
17-16/36-37	0°	0°	0,3 mm	•	6,0 mm	724-082-51	10 pieces
26-27/47-46	0°	0°	0,3 mm	•	6,0 mm	724-083-51	10 pieces
17-16/36-37	0°	6°	0,3 mm	•	6,0 mm	724-086-51	10 pieces
26-27/47-46	0°	6°	0,3 mm	•	6,0 mm	724-087-51	10 pieces
17-16/36-37	0°	6°	0,3 mm	•	4,0 mm	724-092-51	10 pieces
26-27/47-46	0°	6°	0,3 mm	•	4,0 mm	724-093-51	10 pieces
17-16/26-27/47-46/36-37	0°	6°	0,3 mm	•	4,0 mm	724-090-51	10 pieces

CE 0483



Ortho-Cast™ 18 · single rectangular

Slot: 0,46 x 0,64 mm / 18 x 25 · non-convertible · with larger base

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
17-16/36-37	0°	0°	0,3 mm	•	6,0 mm	724-082-50	10 pieces
26-27/47-46	0°	0°	0,3 mm	•	6,0 mm	724-083-50	10 pieces
17-16/36-37	0°	6°	0,3 mm	•	6,0 mm	724-086-50	10 pieces
26-27/47-46	0°	6°	0,3 mm	•	6,0 mm	724-087-50	10 pieces
17-16/36-37	0°	6°	0,3 mm	•	4,0 mm	724-092-50	10 pieces
26-27/47-46	0°	6°	0,3 mm	•	4,0 mm	724-093-50	10 pieces
17-16/26-27/47-46/36-37	0°	6°	0,3 mm	•	4,0 mm	724-090-50	10 pieces

CE 0483

Ortho-Cast™ NF 18 · double rectangular

Slot: 0,46 x 0,64 mm / 18 x 25 · non-convertible · with vertical hook · nickel free

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
47-46/36-37	- 30°	0°	0,3 mm	•	3,5 mm	724-562-51	10 pieces

CE 0483



Ortho-Cast™ 18 · double rectangular

Slot: 0,46 x 0,64 mm / 18 x 25 · non-convertible

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
47-46	- 30°	4°	0,3 mm	•	3,5 mm	724-208-51	10 pieces
36-37	- 30°	4°	0,3 mm	•	3,5 mm	724-209-51	10 pieces

CE 0483



Ortho-Cast™ 18 · double rectangular

Slot: 0,46 x 0,64 mm / 18 x 25 · convertible

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
47-46	- 25°	8°	0,3 mm	•	4,0 mm	724-056-51	10 pieces
36-37	- 25°	8°	0,3 mm	•	4,0 mm	724-057-51	10 pieces
47-46	- 25°	12°	0,3 mm	•	4,0 mm	724-012-51	10 pieces
36-37	- 25°	12°	0,3 mm	•	4,0 mm	724-013-51	10 pieces

CE 0483



Ortho-Cast™ 18 · double rectangular

Slot: 0,46 x 0,64 mm / 18 x 25 · non-convertible · with universal vertical hook

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
17-16/26-27/47-46/36-37	0°	0°	0,3 mm	•	3,5 mm	724-018-51	10 pieces

CE 0483



Ortho-Cast™ 18 · double combination

Slot: 0,46 x 0,64 mm / 18 x 25 · tube 1,15 mm / 45 (occlusal) · non-convertible · for lip bumper / face bow technique

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
17-16	- 14°	14°	0,3 mm	•	4,0 mm	724-204-51	10 pieces
26-27	- 14°	14°	0,3 mm	•	4,0 mm	724-205-51	10 pieces

CE 0483



Ortho-Cast™ NF 18 · triple combination

Slot: 0,46 x 0,64 mm / 18 x 25 · non-convertible · nickel free

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
17-16	- 14°	14°	0,3 mm	•	4,0 mm	724-500-51	10 pieces
26-27	- 14°	14°	0,3 mm	•	4,0 mm	724-501-51	10 pieces

CE 0483



Ortho-Cast™ 18 · triple combination

Slot: 0,46 x 0,64 mm / 18 x 25 · round tube 1,15 mm / 45 (occlusal) · non-convertible · for lip bumper / face bow technique

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
17-16	- 14°	14°	0,3 mm	•	4,0 mm	724-200-51	10 pieces
26-27	- 14°	14°	0,3 mm	•	4,0 mm	724-201-51	10 pieces

CE 0483



Ortho-Cast™ 18 · triple combination

Slot: 0,46 x 0,64 mm / 18 x 25 · round tube 1,15 mm / 45 (occlusal) · convertible · for lip bumper / face bow technique

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
17-16	- 10°	15°	0,3 mm	•	4,0 mm	724-010-51	10 pieces
26-27	- 10°	15°	0,3 mm	•	4,0 mm	724-011-51	10 pieces

CE 0483



For the Bonding Technique



Ortho-Cast™ 22 · single rectangular

Slot: 0,56 x 0,70 mm / 22 x 28 · non-convertible

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
17-16	- 10°	8°	0,3 mm	•	4,0 mm	724-116-51	10 pieces
26-27	- 10°	8°	0,3 mm	•	4,0 mm	724-117-51	10 pieces
17-16	- 14°	10°	0,3 mm	•	4,0 mm	724-214-51	10 pieces
26-27	- 14°	10°	0,3 mm	•	4,0 mm	724-215-51	10 pieces

CE 0483



Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
47-46	- 10°	0°	0,3 mm	•	4,0 mm	724-218-51	10 pieces
36-37	- 10°	0°	0,3 mm	•	4,0 mm	724-219-51	10 pieces
47-46	- 25°	8°	0,3 mm	•	4,0 mm	724-120-51	10 pieces
36-37	- 25°	8°	0,3 mm	•	4,0 mm	724-121-51	10 pieces
47-46	- 30°	12°	0,3 mm	•	4,0 mm	724-114-51	10 pieces
36-37	- 30°	12°	0,3 mm	•	4,0 mm	724-115-51	10 pieces

CE 0483



Ortho-Cast™ 22 · single rectangular

Slot: 0,56 x 0,70 mm / 22 x 28 · convertible

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
47-46	- 20°	0°	0,3 mm	•	4,0 mm	724-228-50	10 pieces
36-37	- 20°	0°	0,3 mm	•	4,0 mm	724-229-50	10 pieces
47-46	- 25°	8°	0,3 mm	•	4,0 mm	724-140-51	10 pieces
36-37	- 25°	8°	0,3 mm	•	4,0 mm	724-141-51	10 pieces

CE 0483



Ortho-Cast™ 22 · single rectangular

Slot: 0,56 x 0,70 mm / 22 x 28 · non-convertible

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
17-16/36-37	0°	0°	0,3 mm	•	6,0 mm	724-084-51	10 pieces
26-27/47-46	0°	0°	0,3 mm	•	6,0 mm	724-085-51	10 pieces
17-16/36-37	0°	6°	0,3 mm	•	6,0 mm	724-088-51	10 pieces
26-27/47-46	0°	6°	0,3 mm	•	6,0 mm	724-089-51	10 pieces
17-16/36-37	0°	6°	0,3 mm	•	4,0 mm	724-094-51	10 pieces
26-27/47-46	0°	6°	0,3 mm	•	4,0 mm	724-095-51	10 pieces
17-16/26-27/47-46/36-37	0°	6°	0,3 mm	•	4,0 mm	724-091-51	10 pieces

CE 0483



Ortho-Cast™ 22 · single rectangular

Slot: 0,56 x 0,70 mm / 22 x 28 · non-convertible · with larger base

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
17-16/36-37	0°	0°	0,3 mm	•	6,0 mm	724-084-50	10 pieces
26-27/47-46	0°	0°	0,3 mm	•	6,0 mm	724-085-50	10 pieces
17-16/36-37	0°	6°	0,3 mm	•	6,0 mm	724-088-50	10 pieces
26-27/47-46	0°	6°	0,3 mm	•	6,0 mm	724-089-50	10 pieces
17-16/36-37	0°	6°	0,3 mm	•	4,0 mm	724-094-50	10 pieces
26-27/47-46	0°	6°	0,3 mm	•	4,0 mm	724-095-50	10 pieces
17-16/26-27/47-46/36-37	0°	6°	0,3 mm	•	4,0 mm	724-091-50	10 pieces

CE 0483



Ortho-Cast™ NF 22 · double rectangular

Slot: 0,56 x 0,70 mm / 22 x 28 · non-convertible · with vertical hook · nickel free

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
47-46/36-37	- 30°	0°	0,3 mm	•	3,5 mm	724-563-51	10 pieces

CE 0483

Ortho-Cast™ 22 · double rectangular

Slot: 0,56 x 0,70 mm / 22 x 28 · non-convertible

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
47-46	- 30°	4°	0,3 mm	•	3,5 mm	724-210-51	10 pieces
36-37	- 30°	4°	0,3 mm	•	3,5 mm	724-211-51	10 pieces



CE 0483

Ortho-Cast™ 22 · double rectangular

Slot: 0,56 x 0,70 mm / 22 x 28 · convertible

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
47-46	- 25°	12°	0,3 mm	•	4,0 mm	724-112-51	10 pieces
36-37	- 25°	12°	0,3 mm	•	4,0 mm	724-113-51	10 pieces



CE 0483

Ortho-Cast™ 22 · double rectangular

Slot: 0,56 x 0,70 mm / 22 x 28 · non-convertible · with vertical hook

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
17-16/26-27/47-46/36-37	0°	0°	0,3 mm	•	3,5 mm	724-019-51	10 pieces



CE 0483

Ortho-Cast™ 22 · double combination

Slot: 0,56 x 0,70 mm / 22 x 28 · tube 1,15 mm / 45 (occlusal) · non-convertible · for lip bumper / face bow technique

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
17-16	- 14°	14°	0,3 mm	•	4,0 mm	724-206-51	10 pieces
26-27	- 14°	14°	0,3 mm	•	4,0 mm	724-207-51	10 pieces



CE 0483

Ortho-Cast™ NF 22 · triple combination

Slot: 0,56 x 0,70 mm / 22 x 28 · non-convertible · nickel free

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
17-16	- 14°	14°	0,3 mm	•	4,0 mm	724-502-51	10 pieces
26-27	- 14°	14°	0,3 mm	•	4,0 mm	724-503-51	10 pieces



CE 0483

Ortho-Cast™ 22 · triple combination

Slot: 0,56 x 0,70 mm / 22 x 28 · round tube 1,15 mm / 45 (occlusal) · non-convertible · for lip bumper / face bow technique

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
17-16	- 14°	14°	0,3 mm	•	4,0 mm	724-202-51	10 pieces
26-27	- 14°	14°	0,3 mm	•	4,0 mm	724-203-51	10 pieces



CE 0483

Ortho-Cast™ 22 · triple combination

Slot: 0,56 x 0,70 mm / 22 x 28 · round tube 1,15 mm / 45 (occlusal) · convertible · for lip bumper / face bow technique

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
17-16	- 10°	15°	0,3 mm	•	4,0 mm	724-110-51	10 pieces
17-16	- 14°	10°	0,3 mm	•	4,0 mm	724-226-50	10 pieces
26-27	- 10°	15°	0,3 mm	•	4,0 mm	724-111-51	10 pieces
26-27	- 14°	10°	0,3 mm	•	4,0 mm	724-227-50	10 pieces



CE 0483

Ortho · single round

round tube 0,9 mm / 36 · non-convertible · for Light Wire Technique · with vertical slot · for lip bumper / face bow technique

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
17-16/36-37	-	6°	0,3 mm	•	6,0 mm	725-033-50	10 pieces
26-27/47-46	-	6°	0,3 mm	•	6,0 mm	725-034-50	10 pieces



CE 0483

For the Bonding Technique



Ortho 18 · single rectangular

Slot: 0,46 x 0,64 mm / 18 x 25 · non-convertible

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
17-16/26-27/47-46/36-37	0°	6°	0,3 mm		4,0 mm	725-310-50	10 pieces

CE 0483



Ortho 18 · single rectangular

Slot: 0,46 x 0,64 mm / 18 x 25 · non-convertible · distal extension

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
17-16/36-37	0°	0°	0,3 mm	•	6,0 mm	725-164-50	10 pieces
26-27/47-46	0°	0°	0,3 mm	•	6,0 mm	725-165-50	10 pieces
17-16/36-37	0°	6°	0,3 mm	•	6,0 mm	725-168-50	10 pieces
26-27/47-46	0°	6°	0,3 mm	•	6,0 mm	725-169-50	10 pieces

CE 0483



Ortho 18 · single rectangular

Slot: 0,46 x 0,64 mm / 18 x 25 · non-convertible

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
17-16/36-37	0°	6°	0,3 mm	•	4,0 mm	725-308-50	10 pieces
26-27/47-46	0°	6°	0,3 mm	•	4,0 mm	725-309-50	10 pieces

CE 0483



Ortho 18 · double combination

Slot: 0,46 x 0,64 mm / 18 x 25 · round tube 1,15 mm / 45 (occlusal) · non-convertible · for lip bumper / face bow technique

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
17-16/26-27/47-46/36-37	0°	0°	0,3 mm		4,6 mm	727-265-50	10 pieces

CE 0483



Ortho 22 · single rectangular

Slot: 0,56 x 0,70 mm / 22 x 28 · non-convertible · distal extension

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
17-16/36-37	0°	0°	0,3 mm	•	6,0 mm	725-166-50	10 pieces
26-27/47-46	0°	0°	0,3 mm	•	6,0 mm	725-167-50	10 pieces
17-16/36-37	0°	6°	0,3 mm	•	6,0 mm	725-170-50	10 pieces
26-27/47-46	0°	6°	0,3 mm	•	6,0 mm	725-171-50	10 pieces

CE 0483



Ortho 22 · single rectangular

Slot: 0,56 x 0,70 mm / 22 x 28 · non-convertible

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
17-16/36-37	0°	6°	0,3 mm	•	4,0 mm	725-306-50	10 pieces
26-27/47-46	0°	6°	0,3 mm	•	4,0 mm	725-307-50	10 pieces
17-16/26-27/47-46/36-37	0°	6°	0,3 mm		4,0 mm	725-311-50	10 pieces

CE 0483



Standard · single round

round tube 0,8 mm / 31

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
17-16/26-27/47-46/36-37	-	0°	0,3 mm		6,0 mm	725-002-50	10 pieces

CE 0483

Standard · single round

round tube 0,9 mm / 36

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
17-16/26-27/47-46/36-37	-	0°	0,3 mm		6,0 mm	725-003-50	10 pieces
17-16/26-27/47-46/36-37	-	0°	0,3 mm		4,6 mm	725-023-50	10 pieces

CE 0483



Standard 18 · double rectangular

Slot: 0,46 x 0,64 mm / 18 x 25 · non-convertible · with vertical hook

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
47-46/36-37	0°	0°	0,3 mm	•	4,6 mm	726-101-50	10 pieces

CE 0483



Rematitan® 18 · double rectangular

Slot: 0,46 x 0,64 mm / 18 x 25 · non-convertible · made of titanium

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
17-16/26-27	- 10°	0°	0,3 mm		4,0 mm	723-020-00	10 pieces

CE 0483



Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
47-46/36-37	- 22°	0°	0,3 mm		4,0 mm	723-022-00	10 pieces

CE 0483



Rematitan® 18 · double rectangular

Slot: 0,46 x 0,64 mm / 18 x 25 · non-convertible · with vertical hook · made of titanium

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
17-16/26-27	- 10°	0°	0,3 mm	•	4,0 mm	723-024-00	10 pieces

CE 0483



Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
47-46/36-37	- 22°	0°	0,3 mm	•	4,0 mm	723-026-00	10 pieces

CE 0483



Rematitan® 22 · double rectangular

Slot: 0,56 x 0,70 mm / 22 x 28 · non-convertible · made of titanium

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
17-16/26-27	- 10°	0°	0,3 mm		4,0 mm	723-021-00	10 pieces

CE 0483



Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
47-46/36-37	- 22°	0°	0,3 mm		4,0 mm	723-023-00	10 pieces

CE 0483



Rematitan® 22 · double rectangular

Slot: 0,56 x 0,70 mm / 22 x 28 · non-convertible · with vertical hook · made of titanium

Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
17-16/26-27	- 10°	0°	0,3 mm	•	4,0 mm	723-025-00	10 pieces

CE 0483



Teeth	Torque	Distal Offset	In/Out	Hook	Width	Order No.	Quantity
47-46/36-37	- 22°	0°	0,3 mm	•	4,0 mm	723-027-00	10 pieces

CE 0483

