



JEWELS

the Funstyle Brackets

NEW · NEW · NEW

The first **colourful** brackets world wide



Trendy colours

Enthusiastic Kids

Youngsters can especially be motivated into orthodontic treatment through bright colours. In removable appliances, it is already a standard option to offer a variety of colours to this target group (for example with Orthocryl® Neon colours and various integrated designs). With JEWELS it is not only possible to offer young patients colourful elastics but also the matching brackets in trendy colours.



Well proven bracket design

Reassurance for both orthodontist and patient

The basis for the new JEWELS brackets are the well established Fascination® ceramic brackets. This means, with JEWELS you benefit from all the construction and material advantages of the Fascination® brackets, including this system's long standing experience. All the colours available in the JEWELS brackets have been tested for bio-compatibility. The result: JEWELS are harmless.



Top Price, Top Product

Optimal value for money

Until now it has only been possible to offer parents metal brackets at a reasonable price for the orthodontic treatment for their children. With JEWELS it is now also possible to conjure up an inexpensive and colourful smile into our children's faces. This has been achieved through the use of first class raw materials, which do not require any specific degree of whiteness and due to the patented production methods used.



JEWELS-Brackets



aquamarin



azur



pink



grenadine

Assortments · Hooks on canines

Roth 18

Roth 22

Teeth	Colour	Amount	REF	REF
13 - 11 / 21 - 23	aquamarin	1 Bracket each = 6 pcs.	700-155-00	700-159-00
13 - 11 / 21 - 23	azur	1 Bracket each = 6 pcs.	700-255-00	700-259-00
13 - 11 / 21 - 23	pink	1 Bracket each = 6 pcs.	700-355-00	700-359-00
13 - 11 / 21 - 23	grenadine	1 Bracket each = 6 pcs.	700-455-00	700-459-00

CE 0483

Single packs · Roth 18

Teeth	Colour	Torque	Angulation	Hooks	Width	REF	Amount
11	aquamarin	+ 11°	+ 5°		3,9 mm	700-115-00	1 piece
21	aquamarin	+ 11°	+ 5°		3,9 mm	700-116-00	1 piece
12	aquamarin	+ 7°	+ 8°		3,4 mm	700-117-00	1 piece
22	aquamarin	+ 7°	+ 8°		3,4 mm	700-118-00	1 piece
13	aquamarin	- 2°	+ 13°	•	3,9 mm	700-121-00	1 piece
23	aquamarin	- 2°	+ 13°	•	3,9 mm	700-122-00	1 piece
11	azur	+ 11°	+ 5°		3,9 mm	700-215-00	1 piece
21	azur	+ 11°	+ 5°		3,9 mm	700-216-00	1 piece
12	azur	+ 7°	+ 8°		3,4 mm	700-217-00	1 piece
22	azur	+ 7°	+ 8°		3,4 mm	700-218-00	1 piece
13	azur	- 2°	+ 13°	•	3,9 mm	700-221-00	1 piece
23	azur	- 2°	+ 13°	•	3,9 mm	700-222-00	1 piece
11	pink	+ 11°	+ 5°		3,9 mm	700-315-00	1 piece
21	pink	+ 11°	+ 5°		3,9 mm	700-316-00	1 piece
12	pink	+ 7°	+ 8°		3,4 mm	700-317-00	1 piece
22	pink	+ 7°	+ 8°		3,4 mm	700-318-00	1 piece
13	pink	- 2°	+ 13°	•	3,9 mm	700-321-00	1 piece
23	pink	- 2°	+ 13°	•	3,9 mm	700-322-00	1 piece
11	grenadine	+ 11°	+ 5°		3,9 mm	700-415-00	1 piece
21	grenadine	+ 11°	+ 5°		3,9 mm	700-416-00	1 piece
12	grenadine	+ 7°	+ 8°		3,4 mm	700-417-00	1 piece
22	grenadine	+ 7°	+ 8°		3,4 mm	700-418-00	1 piece
13	grenadine	- 2°	+ 13°	•	3,9 mm	700-421-00	1 piece
23	grenadine	- 2°	+ 13°	•	3,9 mm	700-422-00	1 piece

CE 0483

More single packs: next page

Single packs · Roth 22

Teeth	Colour	Torque	Angulation	Hooks	Width	REF	Amount
11	aquamarin	+ 11°	+ 5°		3,9 mm	700-135-00	1 piece
21	aquamarin	+ 11°	+ 5°		3,9 mm	700-136-00	1 piece
12	aquamarin	+ 7°	+ 8°		3,4 mm	700-137-00	1 piece
22	aquamarin	+ 7°	+ 8°		3,4 mm	700-138-00	1 piece
13	aquamarin	- 2°	+ 13°	•	3,9 mm	700-141-00	1 piece
23	aquamarin	- 2°	+ 13°	•	3,9 mm	700-142-00	1 piece
11	azur	+ 11°	+ 5°		3,9 mm	700-235-00	1 piece
21	azur	+ 11°	+ 5°		3,9 mm	700-236-00	1 piece
12	azur	+ 7°	+ 8°		3,4 mm	700-237-00	1 piece
22	azur	+ 7°	+ 8°		3,4 mm	700-238-00	1 piece
13	azur	- 2°	+ 13°	•	3,9 mm	700-241-00	1 piece
23	azur	- 2°	+ 13°	•	3,9 mm	700-242-00	1 piece
11	pink	+ 11°	+ 5°		3,9 mm	700-335-00	1 piece
21	pink	+ 11°	+ 5°		3,9 mm	700-336-00	1 piece
12	pink	+ 7°	+ 8°		3,4 mm	700-337-00	1 piece
22	pink	+ 7°	+ 8°		3,4 mm	700-338-00	1 piece
13	pink	- 2°	+ 13°	•	3,9 mm	700-341-00	1 piece
23	pink	- 2°	+ 13°	•	3,9 mm	700-342-00	1 piece
11	grenadine	+ 11°	+ 5°		3,9 mm	700-435-00	1 piece
21	grenadine	+ 11°	+ 5°		3,9 mm	700-436-00	1 piece
12	grenadine	+ 7°	+ 8°		3,4 mm	700-437-00	1 piece
22	grenadine	+ 7°	+ 8°		3,4 mm	700-438-00	1 piece
13	grenadine	- 2°	+ 13°	•	3,9 mm	700-441-00	1 piece
23	grenadine	- 2°	+ 13°	•	3,9 mm	700-442-00	1 piece

CE 0483

Our orthodontic technicians are available to answer any questions you may have regarding our products:

Hotline Tel.-No. Fixed techniques +49 72 31 / 803-550

Hotline Tel.-No. Removable techniques +49 72 31 / 803-555

Fax-No. +49 72 31 / 803-295

Orders by telephone +49 72 31 / 803-210

Further information on Dentaurem products may be found in the internet.

It is possible forthwith to order our products online under:

www.dentaurum.com

Date of information: 01/05


DENTAUREM