

## Premium-Line Pliers / Fixed Appliance Technique



004-350-00	1 piece
------------	---------

CE

### Special Bracket Debonding Pliers, 45°, Angled

#### Sterilizable

For the removal of brackets and bonding adhesive. 45° angled cutting surface. Hard metal inserts. Sharpening possible by the manufacturer.

**Premium-Line · 10 Year Guarantee**



003-349-00	1 piece
------------	---------

CE

### Bracket Removing Pliers

#### Sterilizable

For removing brackets between bracket base and adhesive. Sharpening possible by the manufacturer.

Replacement plastic insert see Accessories Order No. 038-356-00.

**Premium-Line · 10 Year Guarantee**



004-349-00	1 piece
------------	---------

CE

### Bracket Removing Pliers

#### Sterilizable

For removing brackets and adhesive residue. Sharpening possible by the manufacturer. Replacement rubber insert see Accessories Order No. 037-355-00.

**Premium-Line · 10 Year Guarantee**



026-347-00	1 piece
------------	---------

CE

### Oliver-Jones Band Removing Pliers

#### Sterilizable

With replaceable plastic tip. With replaceable hard metal tip. Sharpening possible by the manufacturer. Plastic tip see Accessories Order No. 039-357-00.

**Premium-Line · 10 Year Guarantee**



003-114-00	1 piece
------------	---------



### Johnson Contouring Pliers

#### Sterilizable

With curvature of the beaks only on one side. For contouring and adapting bands and crowns.

**Premium-Line · 10 Year Guarantee**

## Premium-Line Pliers / Fixed Appliance Technique

### Double Beak Pliers


#### Sterilizable

With concave beaks. For forming bicuspid and molar bands in the mouth. Can be used on both sides.

**Premium-Line · 10 Year Guarantee**

015-157-00	1 piece
------------	---------

CE



### Separating Pliers

#### Sterilizable

For expanding Dentalastics® separating rings before positioning.

**Premium-Line · 10 Year Guarantee**

015-160-00	1 piece
------------	---------

CE



### How Pliers, Straight

#### Sterilizable

Universal pliers for the multi-band technique. Due to the slender working tips perfect for placing bands with prewelded brackets, especially in the lower anterior area.

**Premium-Line · 10 Year Guarantee**

003-112-00	1 piece
------------	---------

CE



### How Pliers, Curved

#### Sterilizable

These curved pliers are especially suited for the posterior areas. The curved lip-safe tips adapt especially well to the shape of the tooth.

**Premium-Line · 10 Year Guarantee**

003-111-00	1 piece
------------	---------

CE



### Side Cutter, Maxi

#### Sterilizable

For wires up to  $\varnothing$  1.0 mm / 40, hard. Inserts made of hard metal. Sharpening of inserts possible by the manufacturer.

**Premium-Line · 10 Year Guarantee**

004-266-00	1 piece
------------	---------

CE



### Distal End Cutter, Mini

#### Sterilizable

This distal end cutter has a cut and hold design.

Delicate style. Inserts made of hard metal. Sharpening of inserts possible by the manufacturer.

For wires up to  $\varnothing$  0.56 x 0.70 mm / 22 x 28.

**Premium-Line · 10 Year Guarantee**

**Length: 12 cm / 4 3/4"**

003-701-00	1 piece
------------	---------

CE



## Premium-Line Pliers / Fixed Appliance Technique



003-702-00	1 piece
------------	---------

CE

### Distal End Cutter, Mini

#### Sterilizable

This small size distal end cutter has longer handles plus a cut and hold design. Inserts made of hard metal. Sharpening of inserts possible by the manufacturer.

For wires up to  $\varnothing$  0.56 x 0.70 mm / 22 x 28.

**Premium-Line · 10 Year Guarantee**

**Length: 14 cm / 5 1/2"**



003-700-00	1 piece
------------	---------

CE

### Distal End Cutter, Maxi

#### Sterilizable

This sturdy distal end cutter has a cut and hold design. Inserts made of hard metal. Sharpening of inserts possible by the manufacturer.

For wires up to  $\varnothing$  0.56 x 0.70 mm / 22 x 28.

**Premium-Line · 10 Year Guarantee**

**Length: 12 cm / 4 3/4"**



014-151-00	1 piece
------------	---------

CE

### Ligature Cutter, Mini

#### Sterilizable

Accuracy right to the tip for smooth cutting of soft ligature wires. With hard metal tips replaceable by the manufacturer. Sharpening of inserts possible by the manufacturer.

For wires up to  $\varnothing$  0.4 mm / 16, soft.

**Premium-Line · 10 Year Guarantee**



014-152-00	1 piece
------------	---------

CE

### Ligature Cutter 45°

#### Sterilizable

45° angled cutting tips for easy cutting of ligatures, especially in the posterior areas. With hard metal tips replaceable by the manufacturer. Sharpening of inserts possible by the manufacturer.

For wires up to  $\varnothing$  0.4 mm / 16, soft.

**Premium-Line · 10 Year Guarantee**



016-153-00	1 piece
------------	---------

CE

### Side Cutters, Medium

#### Sterilizable

Black anodized hard metal tips replaceable by the manufacturer. Smooth cuts guaranteed by accurate cutting edges. Sharpening of inserts possible by the manufacturer.

For wires up to  $\varnothing$  0.7 mm / 28, hard.

**Premium-Line · 10 Year Guarantee**

## Premium-Line Pliers / Fixed Appliance Technique

### Light Wire Pliers

#### Sterilizable

For especially fine wire bending work in the light wire technique.

For wires up to  $\varnothing$  0.5 mm / 20, spring hard.

**Premium-Line · 10 Year Guarantee**

003-130-00
------------

1 piece
---------



### Begg Light Wire Pliers

#### Sterilizable

For especially fine wire bending work in the light wire technique. Additional cutting device. Inserts made of hard metal.

For wires up to  $\varnothing$  0.5 mm / 20, spring hard.

**Premium-Line · 10 Year Guarantee**

004-130-00
------------

1 piece
---------



### Angle/Tweed Ribbon Arch Pliers

#### Sterilizable

For applying torque to edgewise wires.

Parallel grip prevents crushing of wire. Hard metal inserts.

For wires up to  $\square$  0.56 x 0.70 mm / 22 x 28.

**Premium-Line · 10 Year Guarantee**

003-442-00
------------

1 piece
---------



### Weingart Universal Pliers

#### Sterilizable

For positioning the arch and bending the arch wire ends. Accurately closing serrated tips. With hard metal inserts replaceable by the manufacturer.

For wires up to  $\varnothing$  0.5 mm / 20, hard.

**Premium Line · 10 Year Guarantee**

CE

003-120-00
------------

1 piece
---------



### Nance Loop Bending Pliers

#### Sterilizable

With four-step beaks for producing loops of different sizes in round and rectangular wires.

For wires up to  $\varnothing$  0.7 mm / 28, hard.

**Premium-Line · 10 Year Guarantee**

000-001-00
------------

1 piece
---------



### Nance Loop Closing Pliers

#### Sterilizable

For bending of stop loops and as special pliers for Adams clasps.

For wires up to  $\varnothing$  0.7 mm / 28, hard.

**Premium-Line · 10 Year Guarantee**

001-002-00
------------

1 piece
---------



## Premium-Line Pliers / Fixed Appliance Technique



003-125-00

1 piece

### Jarabak Light Wire Pliers

#### Sterilizable

For wires up to  $\varnothing$  0.5 mm / 20, spring hard.  
Premium-Line · 10 Year Guarantee



012-350-00

1 piece

### Tweed Loop Forming Pliers

#### Sterilizable

For especially fine bending work.  
For wires up to  $\varnothing$  0.5 mm / 20, hard.  
Premium-Line · 10 Year Guarantee



012-351-00

1 piece

### Tweed Loop Forming Pliers

#### Sterilizable

For especially fine bending work. With individually replaceable hard metal tip. Can be replaced by the user.  
For tips see Accessories  
Order No. 038-357-00.  
For wires up to  $\varnothing$  0.5 mm / 20, hard.  
Premium-Line · 10 Year Guarantee



004-109-00

1 piece

### De La Rosa Contouring Pliers

#### Sterilizable

With guiding grooves for forming round and square arches 16, 18, 22.  
Premium-Line · 10 Year Guarantee



003-185-00

1 piece

### Arch Forming Pliers

#### Sterilizable

For individual construction and activation of lingual/palatal bars. Using these pliers, the bent tips of the wire will fit into all Dentaurem lingual/palatal tubes  $\varnothing$  2 x 0.9 mm / 36.  
For wires up to  $\varnothing$  0.9 mm / 36, hard.  
Premium-Line · 10 Year Guarantee



003-139-00

1 piece

### Angle Wire Bending Pliers, Standard

#### Sterilizable

Sturdy constructed plier with a cone and pyramid shaped beak.  
For wires up to  $\varnothing$  0.7 mm / 28, hard.  
Premium-Line · 10 Year Guarantee

## Premium-Line Pliers / Fixed Appliance Technique

### Angle Wire Bending Pliers

**Sterilizable**

Sturdy constructed plier with a cone and pyramid shaped beak. Inserts made of hard metal.

For wires up to  $\varnothing$  0.7 mm / 28, hard.

**Premium-Line · 10 Year Guarantee**

004-139-00	1 piece
------------	---------



### Loop Forming Pliers, Extra Mini

**Sterilizable**

Ideally suited for forming different loops and adjusting preshaped arches for the fixed appliance technique.

For wires up to  $\varnothing$  0.5 mm / 20, hard.

**Premium-Line · 10 Year Guarantee**

003-350-00	1 piece
------------	---------



### Loop Forming Pliers, Mini

**Sterilizable**

For loop forming and bending of thin round and edgewise wires.

For wires up to  $\varnothing$  0.5 mm / 20, hard.

**Premium-Line · 10 Year Guarantee**

003-351-00	1 piece
------------	---------



### Ligature Forming Pliers

**Sterilizable**

For forming ligatures.

For wires up to  $\varnothing$  0.3 mm / 12, soft.

**Premium-Line · 10 Year Guarantee**

015-158-00	1 piece
------------	---------



### Aderer 3-Prong Pliers, Extra Mini

**Sterilizable**

For bending wires / arches for the fixed appliance technique. Specially suited for bending NiTi wires / arches such as rematitan® "LITE" and Tensic®.

For wires up to max.  $\varnothing$  0.5 mm / 20, hard or  $\varnothing$  0.48 x 0.64 mm / 19 x 25.

**Premium-Line · 10 Year Guarantee**

003-205-00	1 piece
------------	---------

CE



## Premium-Line Pliers / Removable Appliance Technique



003-053-00

1 piece

### Loop Forming Pliers, Medium

#### Sterilizable

For bending small springs and loops.  
For wires up to  $\varnothing$  0.7 mm / 28, hard.

**Premium-Line · 10 Year Guarantee**



003-233-00

1 piece

### Loop Forming Pliers, Maxi

#### Sterilizable

These sturdy pliers have one round and one concave beak for heavier wires, such as space maintainers, etc.

For wires up to  $\varnothing$  0.8 mm / 32, hard.

**Premium-Line · 10 Year Guarantee**



001-650-00

1 piece

### Ruhland Loop Forming Pliers, Magnum

#### Sterilizable

To bend coffin springs, labial arches, loops, etc.

For wires up to  $\varnothing$  1.1 mm / 44, hard.

**Premium-Line · 10 Year Guarantee**



003-515-00

1 piece

### Adams Pliers, Medium

#### Sterilizable

With two smooth, rectangular beaks.  
For wires up to  $\varnothing$  0.7 mm / 28, hard.

**Premium-Line · 10 Year Guarantee**

**Length: 12 cm / 4 3/4"**



004-515-00

1 piece

### Adams Pliers, Medium

#### Sterilizable

With two smooth, rectangular beaks. Inserts made of hard metal.

For wires up to  $\varnothing$  0.7 mm / 28, spring hard.

**Premium-Line · 10 Year Guarantee**

**Length: 12 cm / 4 3/4"**



003-094-00

1 piece

### Adams Pliers, Maxi

#### Sterilizable

With two smooth, rectangular beaks.  
For wires up to  $\varnothing$  0.8 mm / 32, hard.

**Premium-Line · 10 Year Guarantee**

**Length: 14.5 cm / 5 3/4"**

## Premium-Line Pliers / Removable Appliance Technique

### Adams Pliers, Medium

**Sterilizable**

With one smooth, rectangular and one smooth, rounded beak.

For wires up to  $\varnothing$  0.7 mm / 28, hard.

**Premium-Line · 10 Year Guarantee**

**Length: 14 cm / 5 1/2"**

003-516-00	1 piece
------------	---------



### Adams Pliers, Maxi

**Sterilizable**

With one smooth, rectangular and one smooth, rounded beak. Inserts made of hard metal.

For wires up to  $\varnothing$  0.8 mm / 32, hard.

**Premium-Line · 10 Year Guarantee**

**Length: 16 cm / 6 3/4"**

004-095-00	1 piece
------------	---------



### Wire Bending Pliers, Wide

**Sterilizable**

Serrated tips.

Suitable for activating finger springs.

For wires up to  $\varnothing$  0.8 mm / 32, hard.

**Premium-Line · 10 Year Guarantee**

033-452-00	1 piece
------------	---------



### Wire Bending Pliers, Maxi

**Sterilizable**

For bending wires up to  $\varnothing$  0.9 mm / 36, hard.

For cutting wires up to  $\varnothing$  0.7 mm / 28, hard.

With hard metal inserts replaceable by manufacturer.

**Premium-Line · 10 Year Guarantee**

013-522-00	1 piece
------------	---------



### Universal Wire Bending Pliers

**Sterilizable**

With serrated tips and guiding grooves.

For wires up to  $\varnothing$  0.8 mm / 32, hard.

**Premium-Line · 10 Year Guarantee**

003-290-00	1 piece
------------	---------



### Flat Nose Pliers

**Sterilizable**

Smooth faces with guiding grooves. Inserts made of hard metal.

For wires up to  $\varnothing$  0.7 mm / 28, spring hard.

**Premium-Line · 10 Year Guarantee**

004-455-00	1 piece
------------	---------





## Premium-Line Pliers / Removable Appliance Technique



035-455-00

1 piece

### Pointed Beak Pliers

#### Sterilizable

Serrated faces with guiding grooves.  
Fine pointed beaks.

For wires up to  $\varnothing$  0.7 mm / 28, hard.

**Premium-Line · 10 Year Guarantee**



012-074-00

1 piece

### Young Loop Bending Pliers

#### Sterilizable

Universal pliers for different size loops.  
Ideal for labial arches.

For wires up to  $\varnothing$  0.7 mm / 28, hard.

**Premium-Line · 10 Year Guarantee**



003-201-00

1 piece

### Aderer 3-Prong Pliers, Mini

#### Sterilizable

For wires up to  $\varnothing$  0.7 mm / 28, hard.

**Premium-Line · 10 Year Guarantee**

**Length: 12.5 cm / 5"**



003-200-00

1 piece

### Aderer 3-Prong Pliers, Medium

#### Sterilizable

For wires up to  $\varnothing$  0.9 mm / 36, hard.

**Premium-Line · 10 Year Guarantee**

**Length: 13 cm / 5 1/4"**



003-202-00

1 piece

### Aderer 3-Prong Pliers, Maxi

#### Sterilizable

For wires up to  $\varnothing$  1.5 mm / 59, hard.

**Premium-Line · 10 Year Guarantee**

**Length: 15.5 cm / 6 1/4"**



012-520-00


1 piece

### Arrow Clasp Former

#### Sterilizable

For wires up to  $\varnothing$  0.7 mm / 28, hard.

**Premium-Line · 10 Year Guarantee**



## Premium-Line Pliers / Removable Appliance Technique

### Arrow Clasp Bending Pliers


#### Sterilizable

For wires up to  $\varnothing$  0.7 mm / 28, hard.

Premium-Line · 10 Year Guarantee

012-521-00

1 piece



### Face Bow Bending / Cutting Pliers, Magnum

#### Sterilizable

Ideal for bending and cutting the outer bows of a face bow.

For wires up to  $\varnothing$  1.8 mm / 71.

Premium-Line · 10 Year Guarantee

003-180-00

1 piece



## EQ-Line Pliers / Fixed Appliance Technique



045-005-00	1 piece
------------	---------

CE

### Bracket Removing Pliers

#### Sterilizable

For removing brackets.  
Replaceable plastic inserts, see Accessories  
Order No. 038-356-00.



045-021-00	1 piece
------------	---------

CE

### Distal End Cutter

#### Sterilizable

Cut and hold design. Inserts made of hard metal.  
For wires up to  $\varnothing$  0.56 x 0.70 mm / 22 x 28.



045-022-00	1 piece
------------	---------

CE

### Side Cutters, Medium

#### Sterilizable

For cutting wires and arches. Inserts made of hard metal.  
For wires up to  $\varnothing$  0.7 mm / 28, hard.



045-004-00	1 piece
------------	---------

CE

### Oliver-Jones Band Removing Pliers

#### Sterilizable

Replaceable plastic tips.  
Plastic tips see Accessories  
Order No. 039-357-00.



045-006-00	1 piece
------------	---------

CE

### How Pliers, Straight

#### Sterilizable

Universal pliers for the multi-band technique.  
The special tip design makes these pliers suited for seating individual anterior bands.



045-000-00	1 piece
------------	---------

### Angle/Tweed Ribbon Arch Pliers

#### Sterilizable

For applying torque to edgewise wires.  
For wires up to  $\varnothing$  0.56 x 0.70 mm / 22 x 28.

## EQ-Line Pliers / Fixed Appliance Technique

### Weingart Universal Pliers

**Sterilizable**

For positioning the arch and bending the arch wire ends. Accurately closing serrated tips.

For wires up to  $\varnothing$  0.5 mm / 20, hard.

045-003-00	1 piece
------------	---------

CE



### Nance Loop Bending Pliers

**Sterilizable**

With four-step beaks for producing loops of different sizes in round and rectangular wires.

For wires up to  $\varnothing$  0.7 mm / 28, spring hard.

045-010-00	1 piece
------------	---------



### Jarabak Light Wire Pliers

**Sterilizable**

For wires up to  $\varnothing$  0.5 mm / 20, spring hard.

045-009-00	1 piece
------------	---------



### Tweed Loop Forming Pliers

**Sterilizable**

For especially fine bending work.

For wires up to  $\varnothing$  0.5 mm / 20, hard.

045-001-00	1 piece
------------	---------



### Angle Wire Bending Pliers

**Sterilizable**

Sturdy constructed pliers have cone and pyramid shaped beaks.

For wires up to  $\varnothing$  0.7 mm / 28, hard.

045-002-00	1 piece
------------	---------



## EQ-Line Pliers / Removable Appliance Technique



045-012-00

1 piece

### Loop Forming Pliers, Medium

#### Sterilizable

For wires up to  $\varnothing$  0.7 mm / 28, hard.



045-013-00

1 piece

### Loop Forming Pliers, Maxi

#### Sterilizable

These sturdy pliers have one round and one concave beak for heavier wires, such as space maintainers, etc.

For wires up to  $\varnothing$  0.8 mm / 32, hard.



045-014-00

1 piece

### Ruhland Loop Forming Pliers, Magnum

#### Sterilizable

These sturdy pliers have one concave and one round beak for coffin springs, labial arches, etc.

For wires up to  $\varnothing$  1.1 mm / 44, hard.



045-015-00

1 piece

### Adams Pliers, Medium

#### Sterilizable

With two smooth rectangular beaks.

For wires up to  $\varnothing$  0.7 mm / 28, hard.



045-016-00

1 piece

### Adams Pliers, Maxi

#### Sterilizable

With two smooth rectangular beaks.

For wires up to  $\varnothing$  0.8 mm / 32, hard.



045-017-00

1 piece

### Adams Pliers, Maxi

#### Sterilizable

With one smooth, rectangular and one smooth, rounded beak.

For wires up to  $\varnothing$  0.8 mm / 32, hard.

## EQ-Line Pliers / Removable Appliance Technique

### Wire Bending Pliers, Maxi

**Sterilizable**

For bending wires up to  $\varnothing$  0.9 mm / 36, hard. For cutting wires up to  $\varnothing$  0.7 mm / 28, hard.

045-019-00
------------

1 piece
---------



### Flat Nose Pliers

**Sterilizable**

Smooth faces with guiding grooves. Fine pointed beaks. For wires up to  $\varnothing$  0.7 mm / 28, hard.

045-018-00
------------

1 piece
---------



### Young Loop Bending Pliers

**Sterilizable**

Universal pliers for different size loops. Ideal for labial arches. For Wires up to  $\varnothing$  0.7 mm / 28, hard.

045-011-00
------------

1 piece
---------



### Aderer 3-Prong Pliers, Extra Mini

**Sterilizable**

For wires up to  $\varnothing$  0.8 mm / 32, hard.

045-020-00
------------

1 piece
---------



### Aderer 3-Prong Pliers, Mini

**Sterilizable**

For wires up to  $\varnothing$  0.7 mm / 28, hard.

045-007-00
------------

1 piece
---------



### Aderer 3-Prong Pliers, Medium

**Sterilizable**

For wires up to  $\varnothing$  0.9 mm / 36, hard.

045-008-00
------------

1 piece
---------




## Pliers / Removable Appliance Technique



012-523-00

1 piece



### Adams Clasp Bending Pliers

**Can be disinfected but not sterilized.**

Special pliers for easy forming of Adams Clasps in one step.

For wires up to  $\varnothing$  0.7 mm / 28, hard.



037-266-00

1 piece

### Side Cutter, Maxi

**Can be disinfected but not sterilized.**

Made of extra hard vanadium steel, chrome plated. Sharpening of inserts possible by the manufacturer.

For wires up to  $\varnothing$  1.0 mm / 40, hard.



044-177-00

1 piece

### Heavy Gauge Side Cutter

**Can be disinfected but not sterilized.**

Little force required for cutting through large size wires due to special design.

For wires up to  $\varnothing$  1.3 mm / 51, hard.

### Lip Retractor, small size

**Can be disinfected but not sterilized.**  
For efficient work without assistance.  
With flexible plastic connector.

075-601-00	1 piece
------------	---------

CE



### Lip Retractor, standard size

**Can be disinfected but not sterilized.**  
For efficient work without assistance.  
With flexible plastic connector.

075-600-00	1 piece
------------	---------

CE



### Clamping Tweezers

#### Sterilizable

For bracket placement. The blade at the end of the tweezers can be used for angulating edgewise brackets.  
Delicate design.

025-277-00	1 piece
------------	---------

CE



### Clamping Tweezers

#### Sterilizable

With long placement and angulation aid.  
For bracket placement. The blade at the end of the tweezers can be used for angulating edgewise brackets.

025-276-00	1 piece
------------	---------

CE



### Clamping Tweezers

**For posterior teeth.**  
**Sterilizable**

Specially suited for placing brackets and buccal tubes in the posterior region.

025-279-00	1 piece
------------	---------

CE



### Scaler and Band Seater, double ended

#### Sterilizable

Sharp scaling tip on one end for removing adhesive residue.  
Serrated band seater on the other end.

027-349-00	1 piece
------------	---------

CE





## Instruments



023-275-00	1 piece
------------	---------

CE

### Instrument for the Bonding Technique

#### Sterilizable

For smoothing and removing excessive adhesive after positioning of the brackets.



018-000-00	1 piece
------------	---------

CE

### Separating Strip Placer

#### Sterilizable

Holding instrument for separating strips. For smoothing protruding edges of fillings and separating. Complete with 6 separating strips.

Separating strips see Accessories  
Order No. 018-001-00 and 018-002-00.



007-024-00	1 piece
------------	---------

### Scissors, Straight

#### Sterilizable

For cutting thin metal bands.



008-025-00	1 piece
------------	---------

### Scissors, Curved

#### Sterilizable

For cutting thin metal bands.



051-023-00	1 piece
------------	---------

CE

### Eby Band Seater

#### Sterilizable

With curved tip for anterior region.  
Stainless Steel.

Circular snap device.

Seating pressure range 0.8 - 10 N  
(approx. 80 - 1.000 g).

Tips see Accessories Order No. 051-024-00  
and 051-025-00.



026-355-00	1 piece
------------	---------

CE

### Band Seating Instrument, Nylon

#### Sterilizable

With hard metal, triangular serrated tip.  
Can be autoclaved.

### Band Seater

**Sterilizable**  
Large surface tin inlay.

003-006-00	1 piece
------------	---------

CE



### Molar Band Seater

**Sterilizable**  
With replaceable plastic head and hard metal triangular tip.  
Plastic head and hard metal triangular tip see Accessories Order No. 026-357-00.

026-354-00	1 piece
------------	---------

CE



### Mershon Band Pusher

**Sterilizable**  
With serrated rectangular tip and anatomically formed grip for good transfer of force when positioning bands.

024-300-00	1 piece
------------	---------

CE



### Bracket Positioning Gauge

**Sterilizable**  
Simplifies seating the brackets on the teeth with exactly the right distance between bracket slot and incisal edge.  
**With metal tips.**

030-390-00	1 piece
------------	---------

CE



### Bracket Positioning Gauge

**Can be disinfected but not sterilized.**  
Simplifies seating the bracket on the teeth with exactly the right distance between bracket slot and incisal edge. Safety caps for the bracket positioning gauge see Accessories Order No. 030-392-00.  
**With pencils.**

030-391-00	1 piece
------------	---------

CE



### "Mathieu" Needle Holder, Medium

**Sterilizable**  
Ideal for ligating. Convenient and practical design with serrated tips.

000-030-00	1 piece
------------	---------

CE



## Instruments



000-035-00	1 piece
------------	---------

CE

### "Bachmann" Needle Holder with Finger Ring, Medium

#### Sterilizable

Eases ligating with the middle finger under constant pressure. Can also be used in ligating jaw fracture splints. Serrated tips.

**Premium-Line · 10 Year Guarantee**



000-031-00	1 piece
------------	---------

CE

### Needle Holder, Mini

#### Sterilizable

With especially fine serrated tips. For Dentalastics® plastic ligatures.

**Not for wire ligatures.**



000-730-00	1 piece
------------	---------

CE

### Mosquito Forceps, without hook

#### Sterilizable

For placing Dentalastics® plastic ligatures.



000-731-00	1 piece
------------	---------

CE

### Mosquito Forceps, with hook

#### Sterilizable

For placing Dentalastics® plastic ligatures.



023-280-00	1 piece
------------	---------

CE

### Ligator

#### Sterilizable

An instrument for fast ligating of the complete arch. Y-shaped tip simplifies the twisting of the ligatures.

Replacement tips see Accessories

Order No. 023-281-00.

For replacement spools see Accessories

Order No. 023-282-00.

For ligating ligature wires up to  $\varnothing$  0.3 mm / 12, soft.



022-274-00	1 piece
------------	---------

CE

### Ligating Instrument

#### Sterilizable

Ligature adjuster and tucker for long preformed ligatures.

### Ligating Instrument, double ended

#### Sterilizable

For bending and adapting wire ligatures.

021-273-00

1 piece



### Distal End Bender (Langer), double ended

#### Sterilizable

Sturdy instrument for bending arch wire ends.

Handle ø 8 mm

019-000-00

1 piece



### Arch Turret without Torque

Can be disinfected but not sterilized.

For forming flat, plain arches in edgewise wires.

For edgewise wires:

- 0.40 mm / 16
- 0.43 mm / 17
- 0.45 mm / 18
- 0.50 mm / 20
- 0.56 mm / 22

Blue tinted.

004-180-00

1 piece



### Arch Turret with Torque / 16

Can be disinfected but not sterilized.

For forming flat or angled (torque) rectangular arches.

- 0.41 x 0.41 mm / 16 x 16
- in 0°, 7°, 10°, 13° and 16° Torque

Red tinted.

004-186-00

1 piece



### Arch Turret with Torque / 18

Can be disinfected but not sterilized.

For forming flat or angled (torque) rectangular arches.

- 0.46 x 0.64 mm / 18 x 25
- in 0°, 7°, 10°, 13° and 16° Torque

Gold tinted.

004-182-00

1 piece



### Arch Turret with Torque / 22

Can be disinfected but not sterilized.

For forming flat or angled (torque) rectangular arches.

- 0.56 x 0.70 mm / 22 x 28
- in 0°, 7°, 10°, 13° and 16° Torque

Silver tinted.

004-184-00

1 piece



## Instruments



006-013-00	1 piece
------------	---------

CE

### Tension and Compression Gauge

**Can be disinfected but not sterilized.**

Handy device for measuring intraoral tensile and compressive forces. Pressure side with guide slot and fork. Tension side with hook. Measuring scales on pressure and tension sides.

**Measuring range: 28 - 450 g**



040-711-00	1 piece
------------	---------

CE

### Measuring Gauge

**Can be disinfected but not sterilized.**

For measuring tensile and compressive forces.

**Measuring range: 25 - 250 g**



040-712-00	1 piece
------------	---------

CE

**Measuring range: 50 - 500 g**



040-713-00	1 piece
------------	---------

CE

**Measuring range: 100 - 1000 g**

### Dental Vernier (Beerendonk)

**Sterilizable**

With fine, long tips for intraoral measurements.

042-750-00	1 piece
------------	---------



### Dental Vernier (Muenchner Design)

**Sterilizable**

For intra and extraoral measurements, with extra graceful tips.

042-751-00	1 piece
------------	---------



### Orthodontic Ruler (Bernklau)

**Sterilizable**

Made of Stainless Steel. Perfect for all diagnostic measurements, both intraoral and on the model.

044-730-01	1 piece
------------	---------



### Arch Gauge

For measuring labial arches, Adams clasps, palatal bars, etc.

**Flexible white plastic**

766-000-00	1 piece
------------	---------



### Ruler (Muenchner Design)

**Sterilizable**

Stainless Steel, with 1/2 mm and mm increments.

**Length: 15 cm / 6"**

044-731-00	1 piece
------------	---------



### Measuring Grid (Schmuth)

**Can be disinfected but not sterilized.**

This screen is printed red with a 1 x 1 cm scale on one side and the other side in black with a 1 cm x 2 mm scale. When the red and black scales are in alignment, the parallax correction has been achieved for precise measurement of distances.

003-366-00	1 piece
------------	---------



### Orthometer Kit (Korkhaus)

**Can be disinfected but not sterilized.**

Set with orthometer, mm disk, cm disk and space divider.


For space divider see Accessories

Order No. 029-361-01.

029-360-01	1 Set
------------	-------



## Diagnostic



028-353-00

1 piece

CE

### Three Dimensional Bow Divider (Korkhaus)

**Can be disinfected but not sterilized.**

Allows the direction and extent of tooth movement to be checked by a single intra-oral measurement.



030-395-00

1 piece

### Bow Divider

**Can be disinfected but not sterilized.**

Standard design.

To measure tooth width, bracket distance to the incisal edge, etc.




075-400-01

1 Set

### Cephalometric Tracing Kit

**Contents:**

- 4 lead pencils with lead 0.5 mm
- 1 package (= 12 pieces) red refills
- 1 package (= 12 pieces) green refills
- 1 package (= 12 pieces) blue refills
- 1 protractor (ruler)



004-367-00


1 piece

### Diagnostic Ruler-Protractor

For evaluation of radiographs.

**Dimensions: 153 x 38 mm**

**Plastic transparent, clear**



004-368-00


1 piece

### Tracing-Protractor Template

Transparent template with protractor and tooth symbols. For evaluation of radiographs.

**Dimensions: 207 x 132 mm**

**Plastic transparent, clear**



004-369-00

1 piece

### Transparent Template

**Ricketts® method**

For analyzing radiographs. With mm scale, protractor and cut out tooth symbols for tracing.

**Dimensions: 197 x 74 mm**

**Plastic transparent, green**

### Acetate Tracing Foil

For radiograph tracings. One side matte, other side shiny.  
Tear-off pad.

**Dimensions: 18 x 24 cm**

<b>075-100-00</b>	100 pieces
-------------------	------------



**Dimensions: 24 x 30 cm**

<b>075-101-00</b>	100 pieces
-------------------	------------




### Acetate Tracing Foil Cover

For tracing and filing of radiographs.  
One side matte, other side shiny.

**Dimensions: 18 x 24 cm**

<b>075-102-00</b>	50 pieces
-------------------	-----------





## Accessories



055-200-00

50 ml

### Corrosion Protection Oil for Pliers and Instruments

For instrument care before sterilization.



037-355-00

2 pieces

CE

### Rubber Inserts

**Sterilizable**

For Bracket-Removing Plier  
Order No. 004-349-00.



038-356-00

4 pieces

CE

### Plastic Inserts

**Sterilizable**

For Bracket-Removing Pliers  
Order No. 003-349-00 and 045-005-00.



039-357-00

10 pieces

CE

### Plastic Tips

**Sterilizable**

For Oliver-Jones Band Removing Pliers  
Order No. 026-347-00 and 045-004-00.



026-357-00

1 piece

CE

### Plastic Head with triangular hard metal tip

**Sterilizable**

For Molar Band Seater  
Order No. 026-354-00.



038-357-00

3 pieces

### Hard Metal Tip

**Sterilizable**

For Tweed Loop Forming Plier  
Order No. 012-351-00.



030-392-00

4 pieces

CE

### Safety Caps

**Can be disinfected but not sterilized.**

For Bracket Positioning Gauge  
Order No. 030-391-00.

**Color: white**

### Spool

#### Sterilizable

For Ligator Order No. 023-280-00.

023-282-00
------------

1 piece
---------



### Hard Metal Tip

#### Sterilizable

For Ligator Order No. 023-280-00.

023-281-00
------------

1 piece
---------



### Hard Metal Tip, Curved

#### Sterilizable

For Eby Band Seater Order No. 051-023-00.

051-024-00
------------

1 piece
---------



### Hard Metal Tip, Straight

#### Sterilizable

For Eby Band Seater Order No. 051-023-00.

051-025-00
------------

1 piece
---------



### Bow Divider

#### Can be disinfected but not sterilized.

For Korkhaus Orthometer Kit  
Order No. 029-360-01.

029-361-01
------------

1 piece
---------



### Separating Strips, coated on both sides

#### Sterilizable

For Separating Strip Placer  
Order No. 018-000-00.

Width: 4 mm

018-001-00
------------

12 pieces
-----------



## Accessories



018-002-00	12 pieces
------------	-----------

### Separating Strips, coated on one side

#### Sterilizable

For Separating Strip Placer  
Order No. 018-000-00.

Width: 4 mm



055-101-01	1 piece
------------	---------

### Stand for Wires and Pliers, combination

With 6 tubes for straight wires. Do not use disinfectants containing ethanol or alcohol.

Height: Approx. 13 cm / 5 1/8"

Width: Approx. 18 cm / 7"

Made of clear acrylic.



055-000-00	1 piece
------------	---------

### Stand for Pliers

Do not use disinfectants containing ethanol or alcohol.

Length: 13 cm / 5 1/8", for 7 pliers

Made of clear acrylic.



055-100-00	1 piece
------------	---------

Length: 25 cm / 9 7/8", for 14 pliers

Made of clear acrylic.



055-105-00	1 piece
------------	---------

### Stand for Pliers

Do not use disinfectants containing ethanol or alcohol.

For 15 pliers.

Width: 30 cm

Depth: 20 cm

Height: 6.5 cm

Made of white acrylic.

### Bur Stand Magnetic Holder

Made partially of plastic.  
Do not use disinfectants containing ethanol or alcohol.

089-430-00	1 piece
------------	---------



### Bur Stand

Made of plastic.  
Do not use disinfectants containing ethanol or alcohol.

089-431-00	1 piece
------------	---------



### Marking Pencils

#### White Wax Marker

043-700-00	6 pieces
------------	----------



Wipe resistant

#### Black Marker

043-800-00	1 piece
------------	---------



Wipe resistant

#### Red Marker

043-900-00	1 piece
------------	---------

