

TELEREX™



TELEREX™ VIDEO X-RAY FILM VIEWER

- Special CCD & LED technology offers instant viewing of dental X-ray film on a TV or computer monitor.
- View old X-ray film or digitize and archive them to your computer.
- Optional USB Video Grabber enables the Telorex to connect to the computer via the USB port.
- Improves communication between the doctor and patient with enlarged viewing.
- Compact and stylish design with spring loaded film tray, only 9 x 9 x 7 cm.
- Special film adapter included for Pedo size X-ray film.
- 100-240V universal voltage, NTSC or PAL available.

ISO 9001
ISO 13485



DENTAMERICA®

18320 Bedford Circle
P.O. Box 3200
Industry, CA 91744 U.S.A.

Phone: (626) 912-1388
Fax: (626) 913-0510
www.dentamerica.com